



# Skin Pathology Collection Devices

At QML Pathology we are committed to providing the highest quality of service by embracing the latest technology and equipment. In addition to the recent introduction of the punch biopsy device with internal plunger, QML Pathology now has available BIOPBLADES™ for QML Pathology Skin Pathology referrers.

## BIOPBLADE™

The sterile, single-use BIOPBLADE™ is a flexible scalpel used for cutaneous surgery, including: shave biopsy, saucerization of flat lesions and levelling of pedunculated lesions. The unique design of the BIOPBLADE™ incorporates a comfortable and protective 'Fingerguard' in addition to the flexible super sharp blade. This flexibility allows the blade to be positioned at the correct angle for the intended procedure.



The BIOPBLADE™ is utilized for removal of lesions, either elevated (shave biopsy) or flat (saucerization). After the site is anesthetized, the BIOPBLADE™ is held and 'bowed' between the thumb and fingers. The lesion is removed at or just below the surface epithelium. Cosmetic results are normally good and the wound heals without the need for suturing. The Clinician will remove all of the lesion without overly deep penetration to avoid scarring.

## SKIN PUNCH BIOPSY DEVICES

The punch biopsy with internal plunger system allows the lodged skin specimen inside the metal lumen of the punch to be easily ejected.

TYPE OF BIOPSY	SIZES AVAILABLE (MM)
Punch Biopsy with internal plunger	2, 3, 4
Punch Biopsy without plunger	2, 3, 4, 5, 6, 8



Punch Biopsy with internal plunger



Punch Biopsy without plunger

**For further information or to order any of these devices, please contact our Stores Department on (07) 3121 4508 or your local QML Pathology Branch Laboratory.**